


RAYCAST® HIGH PRECISION BASE PLATE-TO-COUCH TOP FIXATION DEVICES

Fixation block for Aktina Couches

Article No. : **32810** (1 pair)



A. GENERAL PRODUCT INFORMATION


The product referred to in these instructions is a device, used to connect High Precision base plates from Orfit Industries to an Aktina Couch with a width of 450mm.

B. PRODUCT DESCRIPTION

This fixation device features two positioning blocks that can be attached on both sides of the long base plates of the AIO and IMRT base plate ranges. Information on these other parts (and instructions on how to make the masks) can be found in the respective 'instructions for use' and on www.orfit.com.

C. POSITIONING ON A BASE PLATE

Mount both blocks onto the side of the base plate using the four M5x25 screws included in the packaging. There are two possible positions on the base plate to mount the blocks.



It is possible to use a set of these fixation blocks together with the Tablefix, or together with a second set of these blocks, to get a better fixation of the base plate onto the Aktina couch. The blocks only fix the base plate in the lateral direction.

D. STORAGE

Always store the product in a safe place to prevent it from getting damaged. Store the product between +10°C (50°F) and 40°C (122°F).

E. PROPERTIES

E.1. Basic Physical Properties

The following physical properties apply to these units:

Dimensions: L 55 mm x W 35 mm x H 42 mm
Weight: 0.054 kg

E.2. Mechanical Properties

This product is made of POM. This provides the necessary stability to get accurate patient positioning even under frequent use.

E.3. Dosimetric Properties

Because of the high density of the POM, do not use this product in the path of the beam.

F. MAINTENANCE AND WASTE MANAGEMENT


This product can be cleaned and disinfected by means of soapy water or an isopropanol based disinfectant, applied with a soft cloth. If unsure about the cleaning fluid, do not use. **Never use aerosol sprays, corrosive cleaning agents, solvents or abrasive detergents.**

Periodic checks of the product should be done to insure the parts are not worn and require repair or replacement. **Do not attempt to make repairs yourself.** Contact your distributor if there are any questions or concerns.

The product can be disposed of with household waste.

G. ADDITIONAL INFORMATION

For additional information such as distributor contact information, product brochures, Safety Data Sheets and regulatory information, please visit our website www.orfit.com.



Note:

It is prohibited to make alterations to this text without prior approval from Orfit Industries. RAYCAST® is a registered trademark of ORFIT Industries.



ORFIT INDUSTRIES
Vosveld 9A | B-2110 Wijnegem | Belgium
T (+32) (0)3 326 20 26
welcome@orfit.com

www.orfit.com



Ref. No. 50000
VERSION 2
LAST UPDATE: 01/04/2021
REVISION DATE: 01/04/2023