

**SAGITTILT™ PRONE BREAST
SOLUTION
SAGITTILT PRONE BREAST
SOLUTION SANITARY COVER**



Article No. : **32070/14**

A. GENERAL PRODUCT INFORMATION

The product referred to in this document is a medical device used for covering the forehead support of the Sagittilt Prone Breast Solution.

B. PRODUCT DESCRIPTION

The Sanitary covers are packaged per 50 pieces in one bag. The Sanitary covers are made of 75% non-woven polypropylene, 20% latex free elastic and 5% thread.

C. USAGE

Place the Sanitary cover over the forehead support of the Sagittilt Prone Breast Solution as shown in the image below. The cover will protect the forehead support and can easily be changed between patients.



D. STORAGE

Always store the product in a safe place to prevent it from getting damaged.

Store the product between +10°C (50°F) and 30°C (86°F).

E. MAINTENANCE AND WASTE MANAGEMENT

This product is designed for single patient use. The product can be disposed of with household waste.

F. ADDITIONAL INFORMATION

For additional information such as distributor contact information, product brochures, Safety Data Sheets and regulatory information, please visit our website www.orfit.com.

Note:

It is prohibited to make alterations to this text without prior approval from Orfit Industries. SAGITTILT™ is a registered trademark of Orfit Industries.



ORFIT INDUSTRIES
Vosveld 9A | B-2110 Wijnegem | Belgium
T (+32) (0)3 326 20 26
welcome@orfit.com

www.orfit.com



Ref. No. 50097
VERSION 2
LAST UPDATE: 02/04/2021
REVISION DATE: 02/04/2023