

## SAGITILT PRONE BREAST SOLUTION SANITARY COVER FOR COMFORT HEAD REST



Article No. : 32070/21

### A. GENERAL PRODUCT INFORMATION

The product referred to in these instructions is a medical device used for covering the comfort head rest of the AIO® Prone Breast Solution and the Sagittilt Prone Breast Solution.

This product may only be used in combination with Orfit immobilisation products. Orfit prohibits the use of unauthorised third-party products in conjunction with its own products.

### B. PRODUCT DESCRIPTION

The Sanitary covers for the comfort head rest are packaged per 50 pieces in one bag. Information on other parts of the AIO® Prone Breast Solution and the Sagittilt Prone Breast Solution can be found in the respective 'instructions for use' and on [www.orfit.com](http://www.orfit.com)

### C. USAGE

Place the Sanitary cover over the comfort head rest as shown in the image below. The cover will protect the comfort head rest and can easily be changed between patients.



### D. PROPERTIES

#### E.1 Physical properties

The following physical properties apply to these sanitary covers :

Dimensions : L 356 mm x W 457 mm

#### E.2 Material properties

Viscose & polyester.

### E. STORAGE

Always store the product in a safe place to prevent it from getting damaged.

Store the product between +10°C (50°F) and 40°C (122°F).

### F. MAINTENANCE AND WASTE MANAGEMENT

This product is designed for single patient use.

The product can be disposed of with household waste.

### G. ADDITIONAL INFORMATION

For additional information such as distributor contact information, product brochures, Safety Data Sheets and regulatory information, please visit our website [www.orfit.com](http://www.orfit.com)

#### Note:

The instructions were written in accordance with the European Directive 93/42/EEC for Medical Devices. It is prohibited to make alterations to this text without prior approval from Orfit Industries. AIO® is a registered trademark of Orfit Industries.



**ORFIT INDUSTRIES**  
Vosveld 9A | B-2110 Wijnegem | Belgium  
T (+32) (0)3 326 20 26  
[welcome@orfit.com](mailto:welcome@orfit.com)

[www.orfit.com](http://www.orfit.com)



Ref. No. 50148  
VERSION 2  
LAST UPDATE: 01/03/2019  
REVISION DATE: 01/03/2021